


# blueGel™

Integrated electrophoresis and visualization system  
QP-1500-01



## TECHNICAL SPECIFICATIONS

### FEATURES

- Gel electrophoresis unit with built-in transilluminator
- Seamless polycarbonate and ABS injection molded parts
- Two-button operation; one for Run, one for Light
- Casting system with integrated comb holder
- Platinum and stainless steel electrodes
- High intensity blue LED transilluminator panel

### GEL SIZE

60 x 60 mm

### GEL VOLUME

15 to 20 ml

### BUFFER VOLUME

30 ml, optimized for TBE

### COMBS

Double sided. For 5, 6, 9 or 13 wells

### POWER SUPPLY

Input: AC 100-240 V, 50-60 HZ  
Output: USB-C PD 2.0+, 12V/1.5A

### DIMENSIONS

9" x 4" x 3" height (23 x 10 x 7 cm.)

### WEIGHT

0.75 lbs., 350 gr

©2020 - minipcr bio, all rights reserved.

Designed in Cambridge, MA | Assembled in the USA | 1770 Massachusetts Ave. Cambridge MA, 02140



+1 781-990-8727



support@minipcr.com



www.minipcr.com

For blueGel™ electrophoresis system orders and buying options visit [minipcr.com](http://minipcr.com)